

# NMR Applications for Polymer Characterisation

B. W. K. Diehl, Spectral Service AG, Emil-Hoffmann-Str. 33, 50996 Cologne, Germany, [diehl@spectralservice.de](mailto:diehl@spectralservice.de)

## The Principles of NMR in Polymer Characterisation


- NMR Chemical Shift instead of Retention Time
- Target molecule must be soluble in deuterated solvent
- Qualitative **and** quantitative analysis is possible.
- NMR signals show molar response.
- One characteristic signal without interference is needed
- Temperature dependent analysis (RT up to 80°C).
- Drugs and/or formulations can be analysed in complex matrices.

Common chromatographical methods often fail in polymer analytics. Therefore multi-nuclear NMR spectroscopy should be established in phar-macopoeia as an alternative for qualitative and quantitative methods.

Polymers are widely spread in drug design.



Exemplary six different target compound groups are presented.

- Silicone (Polydimethylsiloxane)
- Polysaccharides (Drugs and Galenics)
- Polyether (Galenics)
- Polyester (Galenics)
- Poly- and Oligopeptides (Drugs)
- Polyvinyl Compounds (Galenics)





**Fig. 1** shows the  $^1\text{H}$  NMR spectrum of an OH terminated **Polydimethylsiloxane (PDMS)**.

These polymer family can be quantified in drugs using the intensive singlet of the  $\text{Si}(\text{CH}_3)_2$  groups and an internal standard. The method correlates with FT/IR but has a much lower detection limit (aprox. 1ppm). The high spectral dispersion allows the analysis of cyclic types (**Fig. 2**), end-groups (OH, TMS or other modified) and the mean chain length.  $^{29}\text{Si}$  and  $^{13}\text{C}$  satellites are useful calibration signals





**Polysaccharides** represent the biggest group of polymers used in drugs. **Fig. 2** shows the  $^1\text{H}$  NMR spectrum of Aloverose, an acetylated polymannose from Aloe Vera.

The fingerprint signal of the acetate enables an origin test and a quantification. **Fig. 3** shows the comparison of chondroitin sulfates. These high temperature analysis allows to distinguish between different animal origin (e.g. bovine, porcine, shark, chicken). The method is useful for other muco-polysaccharides like heparines, too.



**Polyether, Polyester and Polyvinyl** compounds mostly are used as formulation aids.

**Fig. 4** shows a comparison of  $^{13}\text{C}$  NMR spectra of PEG/PPG copolymers with different sequence of the polymer blocks. **Fig. 5** shows the  $^1\text{H}$  NMR of a drug formulation (top) and for comparison the reference of the corresponding Polymer CoPo Luviscol K 120 (bottom).



**Poly- and oligopeptides** are the most interesting polymers in drugs. Several types already are characterised in some pharmacopoeia. **Fig. 6** shows the  $^1\text{H}$  NMR spectrum of an oligopeptid. The carbonyl region of goserelin is a characteristic part of the  $^{13}\text{C}$  NMR spectrum (see **Fig 7**).

