

# Potential and limitations of quantitative NMR in pharmaceutical analysis

Bernd W.K.Diehl

# Limitations

„Before an analyst uses the method of his college, he will use its toothbrush”

Oder/or

“Was der Bauer nicht kennt, frisst er nicht”

# Potential

NMR gives *molar amounts*

NMR is an *Ab Initio Method*

0,05 % up to 100% Linearity and Dynamic

Selectivity by choosing the nuclei

# Potential

NMR gives *molar amounts*

NMR is an *Ab Initio Method*

0,05 % up to 100% Linearity and Dynamic

Selectivity by choosing the nuclei

NMR is the most powerful detector

Compared to other indirect detection systems

IR, UV, DAD, FID, LSD,

# Methods

- Internal standard
- Twofold analysis using different internal standards
- Standard addition
- External calibration
- External standard

## **SAMPLE PREPARATION!**

- Concentration
- Solvent variation
- Derivatisation
- Temperature

Further upcoming methods (no comments in this presentation):

- HPLC/NMR Coupling
- Solid State NMR


# Applications

- **Definition of standards**  
NMR as a tool to create primary standards. Ab initio determination
- **Identity test for re-import**  
Comparison of active parts and formulation aids
- **Complex mixtures**  
Simultaneous analysis of drugs, by-products, formulation aids and solvents  
Derivatisation, Validation
- **Phytopharmaceuticals**  
Aloe Vera, Lecithin, Fatty oils .....
- **Polymers, Biopolymers**  
Polysaccharide, Chondroitin sulfate, Heparin, modified starch, Peptides,  
Formulation aids, PEG, PVP, silicone ...
- **Heteronuclear NMR**

# Complex Mixtures


## Formulation aids

| Integral            | Initial weight [mg] | Integral | Initial weight [mg] | %          |
|---------------------|---------------------|----------|---------------------|------------|
| TI                  | TI                  | IS       | IS                  | Content    |
| 965,9               | 50,61               | 503,2    | 2,47                | 9,0        |
| 907,9               | 50,08               | 498,3    | 2,60                | 9,1        |
| Molecular weight TI |                     | 397,29   | Mean                | <b>9,1</b> |
| Molecular weight IS |                     | 10       | StdDev              | 0,1        |
| Number of atoms TI  |                     | 5        |                     |            |
| Number of atoms IS  |                     | 197,19   |                     |            |
| Content [%] IS      |                     | 95,7     |                     |            |





## Internal Standard

## Oxaliplatin



# Derivatisation





# Validation


## Reproducibility

|         | Integral            | In. weight [mg] | Integral | In. weight [mg] | mMol  | mMol  | [mg]     | %        | %          |
|---------|---------------------|-----------------|----------|-----------------|-------|-------|----------|----------|------------|
| Test    | P                   | P               | IS       | IS              | IS    | P     | P        | detected | In. weight |
| R1      | 8,7                 | 100,24          | 200,0    | 5,77            | 0,038 | 0,003 | 0,110    | 0,11     | 0,14       |
| R2      | 9,2                 | 99,99           | 200,0    | 5,08            | 0,034 | 0,003 | 0,103    | 0,10     | 0,16       |
| R3      | 9,7                 | 101,09          | 200,0    | 5,81            | 0,039 | 0,004 | 0,124    | 0,12     | 0,15       |
| R4      | 13,1                | 101,28          | 200,0    | 5,53            | 0,037 | 0,005 | 0,159    | 0,16     | 0,15       |
| R5      | 13,4                | 101,22          | 200,0    | 5,91            | 0,039 | 0,005 | 0,174    | 0,17     | 0,17       |
| R6      | 10,6                | 101,15          | 200,0    | 5,94            | 0,040 | 0,004 | 0,138    | 0,14     | 0,15       |
| R7 (L3) | 12,7                | 100,40          | 200,0    | 5,27            | 0,035 | 0,004 | 0,147    | 0,15     | 0,15       |
|         | Molecular weight P  |                 | 33,03    |                 |       |       | Mean     | 0,14     | 0,15       |
|         | Molecular weight IS |                 | 150,22   |                 |       |       | Recovery | 88,1     |            |
|         | Number H (P)        |                 | 3        |                 |       |       | StdDev   | 0,02     |            |
|         | Number H (IS)       |                 | 6        |                 |       |       | StdDev % | 17,2     |            |
|         | Content IS (%)      |                 | 99,9     |                 |       |       |          |          |            |

## Linearity


|      | Integral            | In. weight [mg] | Integral | In. weight [mg] | mMol  | mMol  | [mg]     | %        | %          |
|------|---------------------|-----------------|----------|-----------------|-------|-------|----------|----------|------------|
| Test | P                   | P               | IS       | IS              | IS    | P     | P        | detected | In. weight |
| L1   | 4,1                 | 108,10          | 200,0    | 5,60            | 0,037 | 0,002 | 0,051    | 0,05     | 0,05       |
| L2   | 9,1                 | 101,76          | 200,0    | 5,45            | 0,036 | 0,003 | 0,109    | 0,11     | 0,11       |
| L3   | 12,7                | 100,40          | 200,0    | 5,27            | 0,035 | 0,004 | 0,147    | 0,15     | 0,16       |
| L4   | 23,7                | 100,02          | 200,0    | 5,97            | 0,040 | 0,009 | 0,310    | 0,31     | 0,32       |
| L5   | 47,9                | 105,20          | 200,0    | 5,01            | 0,033 | 0,016 | 0,527    | 0,50     | 0,50       |
| L6   | 85,8                | 107,18          | 200,0    | 5,29            | 0,035 | 0,030 | 0,997    | 0,93     | 0,96       |
| L7   | 165,1               | 101,77          | 200,0    | 5,85            | 0,039 | 0,064 | 2,121    | 2,08     | 2,02       |
|      | Molecular weight P  |                 | 33,03    |                 |       |       |          |          |            |
|      | Molecular weight IS |                 | 150,22   |                 |       |       | Recovery | 98,2     |            |
|      | Number H (P)        |                 | 3        |                 |       |       | StdDev   | 3,6      |            |
|      | Number H (IS)       |                 | 6        |                 |       |       |          |          |            |
|      | Content IS (%)      |                 | 99,9     |                 |       |       |          |          |            |

# Linearity




# Phytopharmaceuticals

Multi Parameter Analysis in Aloe Vera




# Aloverose




# Aloe Vera, Quantification

Internal Standard



# Polymer analysis


PK-Levo



Benserazid = 10.5 %

int. Std.

PVP = 9 %




Integral

200.00

173.91

( ppm )

# PVP, Standard Addition



# Hetero Nuclear NMR

## $^{31}\text{P}$ NMR and $^{19}\text{F}$ NMR

Selectivity, e.g. only phospholipids give signals, triglycerides not

No chromatographic separation needed

**Chemical Shift** instead of RT


$^{13}\text{C}$  Satellites are useful “internal” standards for calibration

## $^{13}\text{C}$ NMR

Analysis of complex mixtures in neat liquid



# $^{31}\text{P}$ -NMR of Soy Lecithin




# <sup>19</sup>F-NMR of Fluazinam

GLB16902-2

- 61.636  
- 61.734  
- 62.081  
- 62.151  
- 62.244


- 63.490  
- 63.560  
- 64.162



<sup>19</sup>F-NMR chemical shifts of by-products

| Chemical shift | Amount [%] |
|----------------|------------|
| -61.63         | 0.04       |
| -61.73         | 0.5        |
| -62.08         | 0.09       |
| -62.15         | 0.09       |
| -62.24         | 0.5        |
| -63.49         | 0.16       |
| -63.56         | 0.13       |
| -64.16         | 0.15       |


## 5-Bach Analysis




- 61.50 - 61.65 - 61.80 - 61.95 - 62.10 - 62.25 - 62.40 - 62.55 - 62.70 - 62.85 - 63.00 - 63.15 - 63.30 - 63.45 - 63.60 - 63.75 - 63.90 - 64.05 - 64.20 - 64.35

( ppm )

# $^{13}\text{C}$ -NMR of Formaldehyde Donor Systems




# $^{13}\text{C}$ -NMR of Formaldehyde



( p p m )


# $^{13}\text{C}$ -NMR of Formaldehyde

97% (w/w) Formaldehyde (50%)




# $^{13}\text{C}$ -NMR of Formaldehyde


97% (w/w) Formaldehyde (55%)



# $^{13}\text{C}$ -NMR of Formaldehyde Donor System



# $^{13}\text{C}$ -NMR of Formaldehyde Donor System





# Conclusion

The potentials of quantitative NMR actually are not conceivable

But, the method must become world wide accepted by the authorities:  
USB, DIN, EU guidelines, Pharmacopoeia etc.

# Conclusion

The potentials of quantitative NMR actually are not conceivable

But, the method must become world wide accepted by the authorities:  
USB, DIN, EU guidelines, Pharmacopoeia etc.

Furthermore you need:

Fantasy, Imagination and Passion

.... and some hundred MHz

# Thanks to:

◆ Spectral Service Team

◆ and to your patience