


## Phospholipid Analysis by $^{31}\text{P}$ -NMR using 600 MHz cryo probe

### Quality Assessment of lecithins 10-fold stronger than common NMR

Quality assessment of soy, milk, egg and marine lecithins is provided by quantitative  $^{31}\text{P}$ -NMR analysis of phospholipids. NMR spectroscopy is the accredited reference method of the I.L.P.S. and was validated in a broad round robin test.

600 MHz  
cQNP  
S/N 4000



### We provide:

Worldwide 48h-Service for Total Quantification of Your Phospholipids!

Individual Reference Standards for HPLC Analysis – created to match your target phospholipid source and composition.

- ✓ native soy lecithin standard SM45
- ✓ milk lecithin standard MK45